

Anti-H_OX40 hlgG2 Antibody(Ivuxolimab)

Product information

GM-23373AB-10	10 µg
GM-23373AB-100	100 µg
GM-23373AB-1000	1 mg

Antibody Information

Species Reactivity	Human
Clone	Ivuxolimab
Source/Isotype	Monoclonal human IgG2, κ
Application	Flow cytometry;
Other Names	ACT35, CD134, IMD16, OX40, TXGP1L
Gene ID	P43489 (human)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

Data Examples

Flow cytometry H_OX40 CHO-K1 Cell Line (Catalog # GM-C18986) was stained with Anti-H_OX40 hlgG2 Antibody(Ivuxolimab) or isotype control antibody.




Fig. FACS

Version:0.3.2.240521