

Anti-TGFB1 hlgG4 Antibody(SRK-181)

Product information


GM-26831AB-10	10 µg
GM-26831AB-100	100 µg
GM-26831AB-1000	1 mg

Antibody Information

Species Reactivity	Human
Clone	SRK-181
Source/Isotype	Monoclonal human IgG4, κ
Application	Flow cytometry;
Other Names	CED, DPD1, IBDIMDE, LAP, TGF-beta1, TGFB, TGFbeta
Gene ID	7040 (human)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

Data Examples

Flow cytometry H_LRRC33 Latent TGF-β1 HEK-293 Cell Line was stained with Anti-TGFB1 hlgG4 Antibody(SRK-181) or isotype control antibody.



SampleID	Geometric Mean : RL1-H
HEK-293 anti-TGFB1+APC-2nd Ab	230
H_LRRC33 Latent TGFB1 HEK-293 H_IgG+APC-2nd Ab	124
H_LRRC33 Latent TGFB1 HEK-293 anti-TGFB1+APC-2nd Ab	6284

Fig. FACS

Version:0.3.2.240521