

Anti-H_CD80 hlgG1 Antibody(Galiximab)

Product information


GM-46075AB-10	10 µg
GM-46075AB-100	100 µg
GM-46075AB-1000	1 mg

Antibody Information

Species Reactivity	Human
Clone	Galiximab
Source/Isotype	Monoclonal human IgG1, λ
Application	Flow cytometry;
Other Names	B7, B7-1, B7.1, BB1, CD28LG, CD28LG1, LAB7
Gene ID	P33681-1 (human)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

Data Examples

Flow cytometry H_CD80 aAPC CHO-K1 Cell Line (Catalog # GM-C24688) was stained with Anti-H_CD80 hlgG1 Antibody(Galiximab) or isotype control antibody.



SampleID	Geometric Mean : FL11-H
CHO-K1 anti-H_CD80+APC-2nd Ab	1916
H_CD80 aAPC CHO-K1 H_IgG+APC-2nd Ab	926
H_CD80 aAPC CHO-K1 anti-H_CD80+APC-2nd Ab	163459

Fig. FACS

Version:0.3.2.240521