

# Anti-TIM3 hlgG2 Antibody(M6903)

## Product information

GM-50627AB-10	10 µg
GM-50627AB-100	100 µg
GM-50627AB-1000	1 mg

## Antibody Information

Species Reactivity	Cynomolgus;
Clone	M6903
Source/Isotype	Monoclonal human IgG2, λ
Application	Flow cytometry;
Other Names	CD366, HAVcr-2, KIM-3, SPTCL, TIMD-3, TIMD3, Tim-3, HAVCR2
Gene ID	A0A2K5U3G7-1 (Cynomolgus)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

## Data Examples

Flow cytometry Cynomolgus\_HAVCR2(TIM3) CHO-K1 Cell Line (Catalog # GM-C25909)was stained with Anti-TIM3 hlgG2 Antibody(M6903) or isotype control antibody.




Fig. FACS

Version:0.3.2.240521