

Anti-MICA/MICB mIgG2a Antibody(7C6)

Product information

GM-51336AB-10	10 µg
GM-51336AB-100	100 µg
GM-51336AB-1000	1 mg

Antibody Information

Species Reactivity	Cynomolgus
Clone	7C6
Source/Isotype	Monoclonal Mouse IgG2a, κ
Application	Flow cytometry;
Other Names	MIC-A; PERB11.1
Gene ID	Q2MGE0-1(cynomolgus)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

Data Examples

Flow cytometry Cynomolgus_MICA(Q2MGE0-1) CHO-K1 Cell Line (Catalog # GM-C26177) was stained with Anti-MICA/MICB mIgG2a Antibody(7C6) or isotype control antibody.




Fig. FACS

Version:0.3.2.240521