

Anti-IL5 hlgG4 Antibody(Reslizumab)

Product information

GM-52408AB-10	10 µg
GM-52408AB-100	100 µg
GM-52408AB-1000	1 mg

Antibody Information

Species Reactivity	Human;
Clone	Reslizumab
Source/Isotype	Monoclonal human IgG4, κ
Application	Block assay;
Other Names	EDF; TRF; IL-5
Gene ID	Q01344-1(human)
Storage	Store at 2-8°C short term (1-2 weeks).Store at ≤ -20°C long term. Avoid repeated freeze-thaw.
Formulation	Phosphate-buffered solution, pH 7.2.
Endotoxin	< 1 EU/mg, determined by LAL gel clotting assay

Data Examples

Block assay

Anti-IL5 hlgG4 Antibody(Reslizumab) inhibits H_IL-5 Reporter 293 Cell Line (Catalog # GM-C06723) Luminescence induced by Interleukin 5 Protein.




Fig. assay