

Biotinylated Cynomolgus NCR3LG1(B7-H6) Protein; His-Avi Tag

Product Information

Product Name Biotinylated Cynomolgus NCR3LG1(B7-H6) Protein; His-Avi Tag

Storage temp. Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.
Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Catalog# / Size **GM-87650RP-25 / 25 μg**
GM-87650RP-200 / 200 μg

Protein Information

Alternative Names B7-H6,NCR3LG1, B7H6,B7 Homolog 6

Source Biotinylated Cynomolgus NCR3LG1(B7-H6) Protein; His-Avi Tag (GM-87650RP) is expressed from human 293 cells (HEK-293). It contains AA Asp 25 - Asp 259 (Accession # XP_005578557.1).

This protein carries a His-Avi tag at the C-terminus.

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/ μg , determined by LAL gel clotting assay

Predicted Mol Mass 28.4 KDa

Formulation Supplied as a 0.2 μm filtered solution of PBS, pH7.4.


Description B7-H6 protein is a co-stimulatory molecule that belongs to the B7 family. It is encoded by the NCR3LG1 gene and is a protein associated with the human immune system. B7-H6 protein was initially discovered in human tumor cells and later detected in other tissues and cells.

B7-H6 protein regulates the activity of NK cells (natural killer cells) by binding to the NKP30 receptor on the surface of NK cells. NK cells are an important type of lymphocyte with the ability to directly kill tumor cells and infected cells, making them crucial members of the immune system.

Research indicates that B7-H6 protein plays a significant role in regulating the activity of NK cells, promoting cytotoxicity, and modulating immune responses. Additionally, the expression of B7-H6 protein is associated with tumor development and immune evasion, making it a potential target for immunotherapy.

Version:3.3 Revision Date:25/12/2023

SDS-PAGE




On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

Bioactivity-ELISA

Bioactivity-ELISA

0.1 μg Biotinylated Cynomolgus NCR3LG1(B7-H6) Protein; His-Avi Tag of per well



Biotinylated Cynomolgus NCR3LG1(B7-H6) Protein; His-Avi Tag (Catalog # GM-87650RP) was immobilized at 1 $\mu\text{g}/\text{ml}$ (100 $\mu\text{L}/\text{well}$) on streptavidin precoated. Increasing concentrations of Anti-B7-H6 hIgG1 Antibody (BI-765049) (Catalog # GM-86750AB) were added.