Genomeditech (Shanghai) Co.,Ltd.
Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

# Mouse MUSK Protein; hFc Tag

#### **Product Information**

Product Name Mouse MUSK Protein; hFc Tag

**Storage temp.** Store at  $\leq$  -70°C, stable for 6 months after receipt.

Recommend to aliquot the protein into smaller quantities for

optimal storage. Please minimize freeze-thaw cycles.

Catalog# / Size GM-87725RP-100 / 100 μg

GM-87725RP-1000 / 1 mg

#### **Protein Information**

Alternative Names Musk, CMS9, FADS, FADS1, muscle associated receptor tyrosine kinase

Source Mouse MUSK Protein; hFc Tag (GM-87725RP) is expressed from human 293

cells (HEK-293). It contains AA Glu22-Thr494 (Accession # Q61006.1).

This protein carries a hFc tag at the C-terminus.

Purity > 95% as determined by SDS-PAGE

Endotoxin < 1 EU/μg, determined by LAL gel clotting assay

Predicted Mol Mass 77.4 KDa

**Formulation** Supplied as a 0.2 μm filtered solution of PBS, pH7.4.

**Description** Musk protein (Muscle-specific Kinase, MUSK) is a tyrosine kinase specifically

expressed in muscle tissue, primarily involved in the development and maintenance of neuromuscular junctions. MUSK plays a key role in synapse formation by promoting the formation and stability of synapses through its

interaction with nerve growth factors released by neurons.

MUSK plays an important role in the development of neuromuscular junctions by interacting with proteins on the neuronal membrane, facilitating the clustering of receptors on the muscle cell membrane. MUSK activates various signaling pathways, regulating the functions and growth of muscle cells. MUSK also plays a crucial role in the differentiation and growth of muscle cells, influencing the

muscle development process.

Abnormalities in MUSK may be associated with certain diseases, such as myasthenia syndrome, making research on MUSK helpful for understanding related physiological and pathological mechanisms. Scientists are exploring the potential of MUSK as a therapeutic target, particularly in the treatment of


neuromuscular diseases.

Genomeditech (Shanghai) Co.,Ltd. Order: 021-68455258/50432826/50432825

Toll-free: 400 627 9288

Email: service@genomeditech.com

### **SDS-PAGE**




On SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

## **Bioactivity-ELISA**

# **Bioactivity-ELISA**

0.2 µg Mouse MUSK Protein; hFc Tag of per well



Mouse MUSK Protein; hFc Tag (Catalog # GM-87725RP) was immobilized at 2  $\mu$ g/ml (100  $\mu$ L/well). Increasing concentrations of Anti-MUSK mIgG2a Antibody(3B2g2m1) (Catalog # GM-87841AB) were added.